



AGFA CR MM3.0 MAMMO CASSETTE 18 X 24

SKU: BA042016 // BB042016 // BC042016 //
BD042016 // BE042016

29,00 €

zzgl. [Versandkosten](#)

Computed radiography (CR) cassettes for diagnostic mammography applications, for use with CR 85-X, CR 35-X digitizers.

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Röntgen & Durchleuchtung](#)

PRODUKTBESCHREIBUNG

Computed radiography (CR) cassettes for diagnostic mammography applications, for use with CR 85-X, CR 35-X digitizers.

Vollfunktionsfähig

Im guten gebrauchten Zustand



Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

