



## [ALOKA UST-939D-3.5 ULTRASOUND TRANSDUCER](#)

**SKU:** 07-032014-01

zzgl. [Versandkosten](#)

Aloka UST-939D-3.5 (3.5 Mhz convex transducer for the  
SSD-650 and 650CL ultrasounds).

Excellent condition

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall / Doppler / Sonden  
/Zubehör](#)

## **PRODUKTBESCHREIBUNG**

Aloka UST-939D-3.5 (3.5 Mhz convex transducer for the SSD-650 and 650CL ultrasounds).

Excellent condition

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



## ZUSÄTZLICHE INFORMATIONEN

