



[DUODENOSCOPE PENTAX ED3401 // ED-3401](#)

SKU: V03-122012

888,00 €

zzgl. [Versandkosten](#)

Vollfunktionsfähig
Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com


info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Flexible Endoskope & Zubehör](#)

PRODUKTBESCHREIBUNG

MODEL **ED-3401**

WORKING CHANNELS1



Diameter, mm (UMBILICAL CORD) 3.8

RECOMMENDED REPROCESSING EtO, glutaraldehyde, peracetic acid

Type (CONFIGURATION) Video duodenoscope

Length, mm (PROBE) 1500

Aspiration (UMBILICAL CORD) Yes

Other functions (Outreach services)

Video signal connector, electrosurgical unit safety cord terminal

EQUIVALENT ILLUMINATION SOURCES/VIDEO PROCESSORS

Pentax EPM-3000/3300 video processor/light source

LENS WASHING (UMBILICAL CORD)

Yes

TIP DEFLECTION RANGE, Å°Up 120, down 90, right 110, left 90

Visual field (OPTICS) 10010° retro

Depth of field, mm (OPTICS) 5

EQUIVALENT CAMERAS (UMBILICAL CORD)

Any video source camera device, video image printer, video slide printer

Length, mm (PROBE) 1250

OUTER DIAMETER, mm (DISINFECTION/STERILIZATION) 11.3

Length markings (INSERTION TUBE)

Every 10 cm

OTHER ATTRIBUTES (Interference compensation) Differenzbesteuerte Ware nach §25a USTG. - **Mwst. wird nicht gesondert ausgewiesen**

High-resolution CCD image chip; 3 remote-control buttons for ancillary devices, VCR, hard copy, image management system; totally immersible; offset biopsy inlet port.

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg



www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN

Weight 10 kg



GALLERIE

