



[EPPENDORF THERMOSTAT PLUS 5352 MIT HEAT BLOCK](#)

SKU: CX112016

799,00 €

zzgl. [Versandkosten](#)

Vollfunktionsfähig
Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com
Mobile 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257



Categories: [Laborequipment](#)

PRODUKTBESCHREIBUNG

This compact, benchtop thermostat maintains different temperature levels up to 99°C, with a cooling feature for temperatures as low as -5°C.

Applications include: Enzyme storage at 0°C, Enzyme reactions at 37°C, Reactions at temperatures typically between 12 and 65°C.

- Temperature and run times are individually programmable
- Stores up to four temperatures and times in succession
- RS-232 interface-provides control and allows documentation in accordance with GLP
- Temperature Control Range of Block from -5 to 99°C
- Range of Possible Temperatures from 30°C below ambient to 99°C
- Minimum Block Temperature: -25°C at ambient of 4°C



- Accuracy of Temperature Control: $\pm 1^{\circ}\text{C}$ between -5°C and 45°C ; $\pm 2^{\circ}\text{C}$ between 45°C and 99°C
- Heating-up Speed: $4^{\circ}\text{C}/\text{min}$
- Cooling-down Speed: $2^{\circ}\text{C}/\text{min}$ between 99°C and 25°C ; $1^{\circ}\text{C}/\text{min}$ between 25°C and -5°C
- Cooling Time: from 22 to 0°C in 20 minutes
- Range of Time Phases: 0:01 to 99:59 hours
- Dimensions, without block (W x D x H): 8.7 x 9.8 x 4.9 in/22 x 25 x 12.5 cm

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN

Weight 12 kg



GALLERIE

