



FUJINON EC-200 MR // EC 200 MR KOLOSKOP

SKU: BI102016-01

zzgl. [Versandkosten](#)

Vollfunktionsfähig
Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com
Mobile 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257



Categories: [Flexible Endoskope & Zubehör](#)

PRODUKTBESCHREIBUNG

Colonoscope

Optical System

- Field of View: 120°
- Direction of View: Forward
- Depth of Field: 5 to100 mm

Insertion Tube

- It Diameter: 13.0 mm

Angulation



-Up/Down: 180°
-Left/Right: 160°

Lenght

-Working Length: 133 cm.
-Total Length: 160 cm.

Biopsy Channel

-Channel Size: 3.2 mm

Compatibility

-Fujinon EPX-201 Video Processor

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN

Weight 10 kg



GALLERIE

