



FUJINON EC-590 WI COLONOSCOPE // KOLONOSKOP

SKU: S022017

4.444,00 €

zzgl. [Versandkosten](#)

Field Of View 140 deg

Depth of Field 3-100 mm

The diameter of the distal end 12.8 mm

The diameter of the flexible part 12.8 mm

Bending Up – 180 deg / down – 180 deg to the
right – 160 deg / left – 160 deg

Working length 1520 mm

Overall Length 1820 mm

Instrument Channel Diameter 3.8mm instrument
channel

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Flexible Endoskope & Zubehör](#)



Field Of View 140 deg
Depth of Field 3-100 mm
The diameter of the distal end 12.8 mm
The diameter of the flexible part 12.8 mm
Bending Up – 180 deg / down – 180 deg to the right – 160 deg / left – 160 deg
Working length 1520 mm
Overall Length 1820 mm
Instrument Channel Diameter 3.8mm instrument channel

Vollfunktionsfähig
Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com

Mobile	0049 176 205 800 86
Tel / Homezone	0049 40 18 20 79 03
Tel / Büro	0049 40 3 61 61 - 256
Fax.	0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN

Weight	11 kg
---------------	-------



GALLERIE

