



**GE ACCESSORIES FOR CATHETE ROOM +  
MULTISYNC LCD 2090UXI MONITOR + HP  
C3980A LASER PRINTER + SYSTEM  
TECHNOLOGIES 545 COMPUTER-SYSTEM**

**SKU:** BI102022

zzgl. [Versandkosten](#)

Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Unkategorisiert](#)

## PRODUKTBESCHREIBUNG

Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

