



**GE CZB 6.5 MHZ PART NR. 2152424 FOR LOGIQ 200**

**SKU: I102014-01**

zzgl. [Versandkosten](#)

CZB Curved Array Transducer 6.5 MHz. For small parts neonatal and biopsy applications.

Part Number 2152424

[Vollfunktionsfähig](#)

[Im guten gebrauchten Zustand](#)

[Artikel-Standort. Hamburg](#)

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall / Doppler / Sonden /Zubehör](#)

## PRODUKTBESCHREIBUNG

CZB Curved Array Transducer 6.5 MHz. For small parts neonatal and biopsy applications.

Part Number 2152424

[Vollfunktionsfähig](#)

[Im guten gebrauchten Zustand](#)



[Artikel-Standort. Hamburg](#)

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile.	0049 176 205 800 86
Tel / Homezone	0049 40 18 20 79 03
Tel / Büro	0049 40 3 61 61 - 256
Fax.	0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

