



KAVO K11 KNEIANLASSER TYP 4980 40000 U/MIN

SKU: 0042015

593,81 €

zzgl. [Versandkosten](#)

KaVo K11 Knieanlasser

Typ 4980

40000 U/min // 200 Watt //

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Dental und Zahnheilkunde](#)

PRODUKTBESCHREIBUNG

KaVo K11 Knieanlasser

Typ 4980

40000 U/min // 200 Watt //

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



ZUSÄTZLICHE INFORMATIONEN



GALLERIE

