



## [MEDISON EC4-9ED ULTRASOUND TRANSDUCER](#)

**SKU:** AK012014-03

zzgl. [Versandkosten](#)

Medison EC4-9/10ED Endo-Cavity Probe  
Specification : 4MHz~9MHz / 150° / 10mm

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)  
0049 40 18 20 79 03  
0049 176 205 800 86  
fax.0049 40 65 86 33 64



---

**Categories:** [Ultraschall / Doppler / Sonden / Zubehör](#)

## **PRODUKTBESCHREIBUNG**

Medison EC4-9/10ED Endo-Cavity Probe  
Specification : 4MHz~9MHz / 150° / 10mm  
Application : OB, Gynecology, Urology

Fits for : Medison SA3200(V1.25.01 or up), SA5000II(V2.12 or up), SA5500, SA6000, SA6000C, SA8800, 128BW, SA8000, PICO, SA9900, X4

Im sehr guten gepflegten Zustand  
Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)  
0049 40 18 20 79 03  
0049 176 205 800 86



fax.0049 40 65 86 33 64

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

