



MEDISON / SONOACE HL5-12ED LINEAR PROBE

SKU: I102015-02

zzgl. [Versandkosten](#)

Medison HL5-12ED Linear Probe

5MHz~12MHz / 40mm / Flat

Application : Vascular, Small part, Breast, Pediatric,
Musculoskeletal

for : Medison Sonoace X6, X8

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall / Doppler / Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

Medison HL5-12ED Linear Probe

5MHz~12MHz / 40mm / Flat

Application : Vascular, Small part, Breast, Pediatric, Musculoskeletal



for : Medison Sonoace X6, X8

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

