



## [MEDTRONIC 7432 NEUROLOGICAL PROGRAMMER WITH CASE](#)

SKU: EG092024

~~299,00 €~~ **234,00 €**

zzgl. [Versandkosten](#)



---

Categories: [Neurologie](#)

## PRODUKTBESCHREIBUNG

The Medtronic 7432 Neurological Programmer is a specialized device used in the field of neuromodulation, particularly for programming and managing implantable neurostimulation devices. It is designed to assist healthcare professionals in adjusting the settings of devices such as deep brain stimulators (DBS) and spinal cord stimulators (SCS).

Key features of the Medtronic 7432 Neurological Programmer may include:

1. **User Interface:** A touchscreen interface that allows for easy navigation and adjustment of stimulation parameters.
2. **Programming Capabilities:** Ability to program various stimulation settings, including frequency, pulse width, and amplitude.
3. **Data Management:** Options for storing patient data and previous programming sessions for reference.
4. **Portability:** Often comes with a protective case for ease of transport and storage.

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

