



[MEDTRONIC SOFAMOR DANEK 906-521](#)
[L032023-K MEDTRONIC SOFAMOR DANEK](#)
[906-521 STRAIGHT PITUITARY RONGEUR](#)
[5MM - 180 - E200 - KK](#)
[STRAIGHT PITUITARY](#)
[RONGEUR 5MM - 180 - E200 - KK](#)

SKU: L032023-K

~~235,00 €~~ **210,60 €**

zzgl. [Versandkosten](#)



Categories: [Unkategorisiert](#)

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

