



**OLYMPUS EXERA TJF-160R VIDEO  
DUODENOSCOPE VIDEODUODENOSCOPE**

SKU: -----

**3.200,00 €**

zzgl. [Versandkosten](#)

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com  
info@mtc-hamburg.com  
Mobile 0049 176 205 800 86  
Tel / Homezone 0049 40 18 20 79 03  
Tel / Büro 0049 40 3 61 61 - 256  
Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Flexible Endoskope & Zubehör](#)

## PRODUKTBESCHREIBUNG

Olympus EXERA TJF-160R Video Duodenoscope.

Removable cap for easy cleaning of the distal end

Superb treatment capability with exceptionally large 4.2mm diameter instrument channel that accommodates a wide range of Endo-Therapy accessories including a variety of stents

Ergonomically designed grip to enhance scope maneuverability and handling

Scope ID function stores individual scope information such as number of procedures performed

Compatible with CV-180, CV-160, CV-140 and CV-100 video processors



Specifications:

Diameter: 11.3mm

Working Length: 124cm

Instrument Channel: 4.2 mm

Field of View: 100 degrees

Angulation Range: Up: 120, Down: 90, Right: 110, Left: 90

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile

0049 176 205 800 86

Tel / Homezone

0049 40 18 20 79 03

Tel / Büro

0049 40 3 61 61 - 256

Fax.

0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

