



OLYMPUS GIF-K10 GASTROSCOPE, GASTROSKOP

SKU: AP092013

zzgl. [Versandkosten](#)

Field Of View : 100°

Insertion Tube : IT Diameter 11.4mm,

Angulation : Up 180°, Down 100°, Right 100°, Left 100°

Length : Working Length 103.5cm

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

0049 40 18 20 79 03

0049 176 205 800 86

fax.0049 40 65 86 33 64



Categories: [Flexible Endoskope & Zubehör](#)

PRODUKTBESCHREIBUNG

Field Of View : 100°

Insertion Tube : IT Diameter 11.4mm,

Angulation : Up 180°, Down 100°, Right 100°, Left 100°

Length : Working Length 103.5cm

Compatibility



Olympus CLV-10 Light Source
Olympus CLV-160 Light Source
Olympus CLV-180 Light Source
Olympus CLV-U20 Light Source
Olympus CLV-U40 Light Source
Olympus CLK-3 Light Source
Olympus CLK-4 Light Source

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

0049 40 18 20 79 03

0049 176 205 800 86

fax.0049 40 65 86 33 64

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

