



OLYMPUS LF-TP ENDOSKOP INTUBATION SCOPE

SKU: AW012016

zzgl. [Versandkosten](#)

- * 1.5mm suction channel
- * 90° field of view
- * 3.8mm distal tip and 4.1mm insertion tube
- * 120° up/down angulation
- * 600mm working length
- * Fully immersible for assured reprocessing
- * Compatible with conventional light sources and video systems for versatility

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Flexible Endoskope & Zubehör](#)

PRODUKTBESCHREIBUNG

- * 1.5mm suction channel



- * 90° field of view
- * 3.8mm distal tip and 4.1mm insertion tube
- * 120° up/down angulation
- * 600mm working length
- * Fully immersible for assured reprocessing
- * Compatible with conventional light sources and video systems for versatility

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

