



OLYMPUS OSF-2 FLEXIBLE FIBER
SIGMOIDOSKOP INK. OLYMPUS LIGHT
SOURCE CABLE LC-2 & KOFFER

SKU: 10-082012-OQ

zzgl. [Versandkosten](#)

gepflegt
Vollfunktionsfähig

www.mtc-hamburg.com
info@mtc-hamburg.com
0049 40 18 20 79 03
0049 176 205 800 86
fax.0049 40 65 86 33 64



Categories: [Flexible Endoskope & Zubehör](#)

PRODUKTBESCHREIBUNG

Olympus OSF-2 Sigmoidoscope -

Flexible Sigmoidoscope with:

- 12.2mm diameter
- 3.2mm channel
- 60cm working length
- 100 Grad field of view
- angulation of 180 Grad /180 Grad (up/down)
- angulation of 160 Grad /160 Grad (right/left).

Differenzbesteuerte Ware nach §25a USTG. - **Mwst. wird nicht gesondert ausgewiesen**



Im gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com
0049 40 18 20 79 03
0049 176 205 800 86
fax.0049 40 65 86 33 64

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

