



[PENTAX EC-3840F2 VIDEO KOLONOSKOP](#)

SKU: BF012016

2.499,00 €

zzgl. [Versandkosten](#)

Field Of View: 120°

IT Diameter : 12.8mm

Angulation Up : 180°

Angulation Down : 180°

Angulation Right : 160°

Angulation Left : 160°

Working Length : 150cm

Channel Size : 3.8mm

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com


info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Flexible Endoskope & Zubehör](#)

PRODUKTBESCHREIBUNG



Field Of View: 120°
IT Diameter : 12.8mm
Angulation Up : 180°
Angulation Down : 180°
Angulation Right : 160°
Angulation Left : 160°
Working Length : 150cm
Channel Size : 3.8mm

Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com

Mobile. 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

