



## [PENTAX ED-3430K DUODENOSCOPE / DUDONOSKOP](#)

**SKU:** AI092013 - AJ092013

**2.499,00 €**

zzgl. [Versandkosten](#)

Field Of View 120°

Direction Of View : Side

Depth Of Field : 4 - 100 mm

Insertion Tube : IT Diameter 11.3mm ,

Angulation : Up 120°, Down 90 °, Right 110 °, Left 90 °

Length: Working Length 125cm, Total Length 157cm

Biopsy Channel: Channel Size 3.8 mm

Compatibility : EPK-700, EPK-1000

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Flexible Endoskope & Zubehör](#)

## PRODUKTBESCHREIBUNG

Field Of View 120°



Direction Of View : Side

Depth Of Field : 4 - 100 mm

Insertion Tube : IT Diameter 11.3mm ,

Angulation : Up 120°, Down 90°, Right 110°, Left 90°

Length: Working Length 125cm, Total Length 157cm

Biopsy Channel: Channel Size 3.8 mm

Compatibility : EPK-700, EPK-1000

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile

0049 176 205 800 86

Tel / Homezone

0049 40 18 20 79 03

Tel / Büro

0049 40 3 61 61 - 256

Fax.

0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN

**Weight**

18 kg



## GALLERIE

