



PENTAX EG-3430K GASTROSKOP / GASTROSCOPE

SKU: AM092013

zzgl. [Versandkosten](#)

Field Of View : 140°

Direction Of View : Forward

Depth Of Field : 4-100mm

Insertion Tube : IT Diameter 11.4mm,

Angulation : Up 210°, Down 120 °, Right 120 °, Left 120 °

Length : Working Length 105 cm, Total Length 137.7cm

Biopsy Channel : Channel Size 3.8mm

gepflegt

Inkl.Koffer

Vollfunktionsfähig

www.mtc-hamburg.com

info@mtc-hamburg.com

0049 40 18 20 79 03

0049 176 205 800 86

Fax.0049 40 65 86 33 64



Categories: [Flexible Endoskope & Zubehör](#)



PRODUKTBESCHREIBUNG

Field Of View : 140°

Direction Of View : Forward

Depth Of Field : 4-100mm

Insertion Tube : IT Diameter 11.4mm,

Angulation : Up 210°, Down 120 °, Right 120 °, Left 120 °

Length : Working Length 105 cm, Total Length 137.7cm

Biopsy Channel : Channel Size 3.8mm

gepflegt

Inkl.Koffer

Vollfunktionsfähig

www.mtc-hamburg.com

info@mtc-hamburg.com

0049 40 18 20 79 03

0049 176 205 800 86

Fax.0049 40 65 86 33 64

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

