



SHANDON CONSUL AUTOMATED ROBOTIC COVERSLIPPER

SKU: ER042020-V

2.378,81 €

zzgl. [Versandkosten](#)

Consul is a robotic, fully automated glass coverslipping instrument capable of up to 400 slides per hour."

220v, 50/60hz,

Zustand. Gebrauchte

Artike-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Labor](#), [Laborequipment](#)

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

