



## SIEMENS 5MHZ LINEAR SONDE / PROBE

**SKU:** J032013 & 03-062014 // AJ042014-02

zzgl. [Versandkosten](#)

für Sonoline SX/SL/SL1/SL2/CD AC/Si 200/Si250 Si 400 /Si 450

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall / Doppler / Sonden](#)  
[/Zubehör](#)

## PRODUKTBESCHREIBUNG

für Sonoline SX/SL/SL1/SL2/CD AC/Si 200/Si250 Si 400 /Si 450

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

