



## SONOACE L5-9 PB-L75/60ED LINEAR ARRAY ULTRASOUND TRANSDUCER PROBE

**SKU:** F032014-02

zzgl. [Versandkosten](#)

geeignet, unter anderen für .SONOACE // MEDISON  
SA-5500; 6000C; 6000-II ; 8000 ...

Frequency 5.0-9 Mhz

Anwendung. Abdomen, Fetal heart, Gynecology and OB

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall / Doppler / Sonden](#)  
[/Zubehör](#)

## **PRODUKTBESCHREIBUNG**

geeignet, unter anderen für .SONOACE // MEDISON SA-5500; 6000C; 6000-II ; 8000 ...

Frequency 5.0-9 Mhz

Anwendung. Abdomen, Fetal heart, Gynecology and OB

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86



Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

