



[SYSMEX AD-270 // AD 270 AUTO DILUTOR](#)

SKU: I042016-VER

zzgl. [Versandkosten](#)

Sysmex Auto Dilutor Type AD-270 is used for diluting blood samples in 2 stages:

1 stage - for WBC determination

2 stage - for RBC determination

Blood is diluted in dilutor with isotonic solution and next WBC, RBC HCT, HBG and derived parameters are measured in Cell Counter.

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Laborequipment](#)

PRODUKTBESCHREIBUNG

Sysmex Auto Dilutor Type AD-270 is used for diluting blood samples in 2 stages:

1 stage - for WBC determination

2 stage - for RBC determination

Blood is diluted in dilutor with isotonic solution and next WBC, RBC HCT, HBG and derived parameters are measured in Cell Counter.

Vollfunktionsfähig

Im guten gebrauchten Zustand



Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN

Weight

7 kg



GALLERIE

