



## [SYSMEX AUTO DILUTOR AD-260](#)

SKU: V042016

999,00 €

zzgl. [Versandkosten](#)

Sysmex AD-260 – Auto Dilutor

Manufactures: Toa, Medical Electronics Co. Ltd.,  
Kobe, Japan

Sysmex Auto Dilutor Type AD-260 is used for  
diluting blood samples in 2 stages:

1 stage – for WBC determination

2 stage – for RBC determination

Blood is diluted in dilutor with isotonic solution and  
next WBC, RBC HCT, HBG and derived parameters  
are measured in Cell Counter.

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

Categories: [Laborequipment](#)

## PRODUKTBESCHREIBUNG

Sysmex AD-260 – Auto Dilutor

Manufactures: Toa, Medical Electronics Co. Ltd., Kobe, Japan



Sysmex Auto Dilutor Type AD-260 is used for diluting blood samples in 2 stages:

1 stage – for WBC determination

2 stage – for RBC determination

Blood is diluted in dilutor with isotonic solution and next WBC, RBC HCT, HBG and derived parameters are measured in Cell Counter.

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN

**Weight**

5 kg



## GALLERIE

