



**TOSHIBA PVG-366M / PVG 366 M 3,75MHZ,  
KONVEX PROBE**

**SKU:** BW022016 // AS052018

**799,00 €**

zzgl. [Versandkosten](#)

Toshiba Konvex 3,75 MHz

PVG-366M

Konvex 3.75MHZ

Anwendung : Abdomen

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall / Doppler / Sonden](#)  
[/Zubehör](#)

## PRODUKTBESCHREIBUNG

Toshiba Konvex 3,75 MHz

PVG-366M

Konvex 3.75MHZ

Anwendung : Abdomen

Artikel-Standort. Hamburg



[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile.            0049 176 205 800 86  
Tel / Homezone    0049 40 18 20 79 03  
Tel / Büro         0049 40 3 61 61 - 256  
Fax.                0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

