



## WELCH ALLYN PIC40 E BIPHASIC AED TRAINING DEFIBRILLATOR

**SKU:** Q112016 // AD112016

**888,00 €**

zzgl. [Versandkosten](#)

Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg

www.mtc-hamburg.com  
info@mtc-hamburg.com  
Mobile 0049 176 205 800 86  
Tel / Homezone 0049 40 18 20 79 03  
Tel / Büro 0049 40 3 61 61 - 256  
Fax. 0049 40 3 61 61 - 257



---

Categories: [Defibrillatoren](#)

## PRODUKTBESCHREIBUNG

COLOR SCREEN

EXTERNAL PADDLES SET AND PEDIATRIC ADAPTER

Medicine calculations

NON INVASIVE BLOOD PRESURE WITH CUFF AND HOSE

CO2 ( CAPNOGRAPH FEATURE) OPTION INCLUDED NO CABLE OR TUBING INCLUDED

OXYGEN SATURATION WITH ORIGINAL NELLCOR SENSOR AND WELCH ALLYN EXTENSION CABLE

ECG WITH 3 LEADS CABLE

PACEMAKER/ AED WITH ADAPTER

PRINTER WITH PAPER

DEFIBRILATION/EXTERNAL PACING VERIFICATION DEVICE

BATTERY AND POWER ADAPTER INCLUDED (100-220 volts)



Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN

**Weight** 15 kg



## GALLERIE

